



MEDICAL INDUSTRIES

Hoffmann Family of Companies



NOA Platinum SC Plus Bed

Quality Features at an Affordable Price

Features

- Patented in bed scale system
- Scale never needs calibration
- Patient hand pendant
- Transport handles for ease of transporting bed and patient
- Digital nurse panel providing you real time readout of head deck angle, knee deck angle, height of bed and angle of T/RT
- 15 year warranty on deck and frame welds
- 4 year warranty on electronics
- One button CPR quick release head deck
- Central and directional locking system with 5" casters
- Bed exit alarm with multiple zones that ties into the nurse call system
- Underbed lighting
- Optional side rail controls for clinician ease of use
- 600 lb. safe working load
- Meets 60601-2-52



Technical Specifications

NOA Platinum SC Plus Hospital Bed

Frame width	35"	Safe working load	600 lbs.
Frame length with head and footboard		Deck height range	11" - 31"
Normal length	84.5"	Max & reverse trendelenburg	14°
Sleep deck		Caster size	5"
Width	35"	Max back angle	70°
Length	80"	Max hip-knee angle	35°
Bed clearance	5.5"		
System			
Accuracy	+/- 1% of patient weight		
Weigh in any position	Yes		
Scale repeatability	2% of patient weight		

